

PureLink™ RNA Mini Kit

Catalog No. 12183020

Quantity: 10 preps

Store at room temperature

Contents and Storage

The components of the PureLink™ RNA Mini Kit are shipped at room temperature. Upon receipt, store all components at room temperature. The components included with the kit are listed in the table below. Sufficient reagents are provided to perform 10 preparations.

Components	Quantity
Lysis Buffer	20 mL
Wash Buffer I	10 mL
Wash Buffer II	4 mL
RNase-Free Water	3 mL
Spin Cartridges (with collection tubes)	10 each
Collection Tubes	10 each
Recovery Tubes	10 each

Description

The PureLink™ RNA Mini Kit provides a simple, reliable, and rapid method for isolating high-quality total RNA from a wide variety of samples, including animal and plant cells and tissue, blood, bacteria, yeast, and liquid samples. The purified total RNA is suitable for use in a variety of downstream applications (see page 4).

Note: This manual provides protocols for purifying RNA from animal and plant cells only. The RNA Mini Kit allows you to purify RNA from a variety of samples and using alternate protocols. **If you wish to purify RNA from other samples or use these alternate protocols, refer to the main RNA Mini Kit manual**, available for downloading from the product page web site at www.invitrogen.com or by contacting Technical Support.

System Overview

In general, samples are lysed and then homogenized in the presence of guanidinium isothiocyanate, a chaotropic salt capable of protecting the RNA from endogenous RNases (Chirgwin *et al.*, 1979). After homogenization, ethanol is added to the sample. The sample is then processed through a Spin Cartridge containing a clear silica-based membrane to which the RNA binds. Any impurities are effectively removed by subsequent washing (Vogelstein & Gillespie, 1979). The purified total RNA is then eluted in RNase-Free Water (Tris Buffer, pH 7.5 may also be used).

Kit Specifications

Cartridge Binding Capacity:	~1 mg nucleic acid
Cartridge Reservoir Capacity:	700 µL
Wash Tube Capacity:	2.0 mL
Centrifuge Compatibility:	Capable of centrifuging >12,000 × g
Elution Volume:	30 µL–3 × 100 µL (3 sequential elutions with 100 µL each)
RNA Yield	Varies with sample type and quality

General Guidelines

- The maximum RNA binding capacity of the PureLink™ RNA Mini Kit Spin Cartridge is ~1 mg. If you are processing samples that contain more than 1 mg of total RNA, divide the sample into aliquots such that each aliquot contains <1 mg total RNA for each Spin Cartridge used.
- For an optional protocol with on-column DNase treatment, refer to the full PureLink™ RNA Mini Kit manual.
- Use disposable, individually wrapped, sterile plastic ware and sterile, disposable RNase-free pipette tips and tubes.
- Wear disposable gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin; change gloves frequently, particularly as the protocol progresses from crude extracts to more purified material.
- Always use proper microbiological aseptic techniques when working with RNA.
- When purifying total RNA from fresh samples, keep fresh cell and tissue samples on ice immediately after harvesting and quickly proceed to **Lysis and Homogenization** (page 2).
- When purifying total RNA from frozen samples, freeze samples immediately after collection in liquid nitrogen or on dry ice. Keep samples frozen at –80°C or in liquid nitrogen until proceeding to **Lysis and Homogenization** (page 2).
- Use RNase AWAY® Reagent (cat. no. 10328-011) to remove RNase contamination from work surfaces and non-disposable items such as centrifuges and pipettes used during purification.

Caution

- Both Lysis Buffer and Wash Buffer I contain guanidine isothiocyanate (an irritant). This chemical is harmful when in contact with the skin, or when it is inhaled or ingested. **Do not** add bleach or acidic solutions directly to solutions or sample preparation waste that contains guanidinium isothiocyanate, as reactive compounds and toxic gases are formed.
- Solutions containing ethanol are considered flammable. Use appropriate precautions when using this chemical.
- For your protection, always wear a laboratory coat, gloves and safety glasses when handling these chemicals. Dispose of the buffers and chemicals in appropriate waste containers.

Purifying RNA from Animal and Plant Cells

Materials Needed

You will need the following items in addition to the kit components:

- 2-mercaptoethanol
- 70% ethanol (in RNase-Free Water)
- 1.5 mL RNase-free microcentrifuge tubes
- PBS ($>10^7$ cells per sample)
- Homogenizer, or RNase-free syringe (1 mL) with 18–21-gauge needle or, Rotor-stator homogenizer
- 96–100% ethanol
- Microcentrifuge capable of centrifuging $12,000 \times g$
- 15 mL RNase-free tubes ($>10^7$ cells per sample)
- RNase-free pipet tips

Buffer Preparation and Parameters

Preparing Wash Buffer II with Ethanol

Before beginning lysis and homogenization, add 16 mL 96–100% ethanol to Wash Buffer II. Check the box on the Wash Buffer II label to indicate that ethanol was added. Store Wash Buffer II with ethanol at room temperature.

Preparing Lysis Buffer with 2-Mercaptoethanol

Before beginning lysis and homogenization, prepare a fresh amount of Lysis Buffer containing 1% 2-mercaptoethanol for each purification procedure. *Add 10 μ L 2-mercaptoethanol for each 1 mL Lysis Buffer used.*

Amount of Lysis Buffer Needed

Number of cells in your sample	Amount of Lysis Buffer Needed (prepared with 2-mercaptoethanol)
$\leq 1 \times 10^6$	0.3 mL* (*Use 0.6 mL if using rotor-stator for lysis or homogenization)
1×10^6 – 5×10^6	0.6 mL
5×10^6 – 5×10^7	0.6 mL per 5×10^6 cells (e.g., use 1.2 mL for 1×10^7 cells and 6.0 mL for 5×10^7 cells)

Lysis and Homogenization

$\leq 5 \times 10^6$ Suspension Cells

- Transfer cells to an RNase-free tube and centrifuge at $2,000 \times g$ for 5 minutes at 4°C to pellet. Discard growth medium.
- Add 0.3 **or** 0.6 mL Lysis Buffer with 2-mercaptoethanol to your sample (see table above for correct amounts).
- Vortex at high speed until the cell pellet is completely dispersed and the cells appear lysed.
- Proceed with **one** of the following homogenization options at room temperature:
 - Transfer the lysate to a Homogenizer inserted in a Collection Tube and centrifuge at $12,000 \times g$ for 2 minutes, **or**
 - Pass the lysate 5–10 times through an 18- to 21-gauge needle attached to an RNase-free syringe, **or**
 - Transfer the lysate to an RNase-free tube and homogenize using a rotor-stator homogenizer at maximum speed for at least 45 seconds. Centrifuge the homogenate at $\sim 2,600 \times g$ for 5 minutes, then transfer the supernatant to a clean RNase-free tube.
- Proceed to **Binding, Washing, and Elution**, next page.

$\leq 5 \times 10^6$ Monolayer Cells

- Remove the growth medium from the cells, then add 0.3 or 0.6 mL Lysis Buffer with 2-mercaptoethanol to your sample (see table above for correct amounts).
- Proceed with **one** of the following homogenization options at room temperature:
 - Transfer the lysate to a Homogenizer inserted in a Collection Tube and centrifuge at $12,000 \times g$ for 2 minutes, **or**
 - Transfer the lysate to a 1.5 mL tube and pass 5–10 times through an 18- to 21-gauge needle attached to a syringe, **or**
 - Transfer the lysate to an appropriately RNase-free tube and homogenize using a rotor-stator homogenizer at maximum speed for at least 45 seconds. Centrifuge the homogenate at $\sim 2,600 \times g$ for 5 minutes, then transfer the supernatant to a clean RNase-free tube.
- Proceed to **Binding, Washing, and Elution**, next page.

Lysis and Homogenization, Continued

5×10^6 – 5×10^7 Suspension Cells

Note: This protocol uses a rotor-stator homogenizer. Use 15-mL RNase-free tubes to compensate for volume expansion.

1. Transfer cells to a 15-mL tube and centrifuge at $2,000 \times g$ for 5 minutes at 4°C to pellet. Remove and discard the supernatant.
2. Add 0.6 mL Lysis Buffer with 2-mercaptoethanol per 5×10^6 cells to your sample (refer to previous page for amounts).
3. Vortex at high speed until the cell pellet is completely dispersed and the cells appear lysed.
4. Homogenize cells using a rotor-stator homogenizer at maximum speed for at least 45 seconds, then centrifuge the homogenate at $\sim 2,600 \times g$ for 5 minutes at room temperature then transfer the supernatant to a clean 15-mL RNase-free tube.
5. Proceed to **Binding, Washing, and Elution**, below

Frozen Cell Pellets

1. Transfer the frozen cell pellet to an RNase-free tube then **add** the appropriate volume of Lysis Buffer prepared with 2-mercaptoethanol to your sample (refer to the previous page for correct amounts).
2. Vortex at high speed until the cell pellet is completely dispersed and the cells appear lysed.
3. Proceed with **one** of the following homogenization options at room temperature:
 - a. Transfer the lysate to a Homogenizer inserted in a Collection Tube and centrifuge at $12,000 \times g$ for 2 minutes, **or**
 - b. Pass the lysate 5–10 times through an 18- to 21-gauge needle attached to an RNase-free syringe, **or**
 - c. Transfer the lysate to an RNase-free tube and homogenize using a rotor-stator homogenizer at maximum speed for at least 45 seconds. Centrifuge the homogenate at $\sim 2,600 \times g$ for 5 minutes, then transfer the supernatant to a clean RNase-free tube.
4. Proceed to **Binding, Washing, and Elution**, below.

Binding, Washing and Elution

Follow the steps below to bind, wash, and elute the RNA from your sample:

1. Add one volume 70% ethanol to each volume of cell homogenate (prepared as described above).
Note: If part of the sample was lost during homogenization, adjust the volume of ethanol accordingly.
2. Vortex to mix thoroughly and to disperse any visible precipitate that may form after adding ethanol.
3. Transfer up to 700 μL of the sample (including any remaining precipitate) to the Spin Cartridge (with the Collection Tube).
4. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through, and reinsert the Spin Cartridge into the same Collection Tube.
Note: When extracting RNA from larger quantities of starting material the homogenate may become viscous, requiring additional centrifugation (up to 15 minutes) for the sample to pass completely through the column.
5. **Repeat** Steps 3–4 until the entire sample has been processed.
6. Add 700 μL Wash Buffer I to the Spin Cartridge.
7. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through **and** the Collection Tube. Place the Spin Cartridge into a **new** Collection Tube.
8. Add 500 μL Wash Buffer II with ethanol (previous page) to the Spin Cartridge.
9. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through.
10. Repeat Steps 8–9 **once**.
11. Centrifuge the Spin Cartridge at $12,000 \times g$ for 1–2 minutes to dry the membrane with attached the RNA. **Discard** the Collection Tube and insert the Spin Cartridge into a Recovery Tube.
12. Add 30 μL –100 μL RNase-Free Water to the center of the Spin Cartridge.
13. Incubate at room temperature for 1 minute.
14. Centrifuge the Spin Cartridge for 2 minutes at $\geq 12,000 \times g$ at room temperature to elute the RNA from the membrane into the Recovery tube.
Note: If the expected RNA yield is $>100 \mu\text{g}$, perform 3 sequential elutions of 100 μL each. Collect the eluates in a single tube.
15. Store your purified RNA or proceed to downstream application of choice.

Storage and Downstream Applications of Purified RNA

Store your purified RNA on ice if you will use the RNA within a few hours. For long-term storage, store your purified RNA at -80°C .

You may use the purified total RNA for qRT-PCR, northern blotting, nuclease protection assays, RNA amplification for microarray analysis, cDNA library preparation after poly(A)⁺ selection, or any desired downstream application.

If highly pure RNA without genomic DNA contamination is required, perform DNase I treatment after purification (refer to the PureLink™ RNA Mini Kit manual available from our web site at www.invitrogen.com for details). You can determine the quality and quantity of the purified RNA using UV absorbance at 260 nm or with the Quant-iT™ RNA Assay Kit (Cat. no. Q33140).

Accessory Products

A large selection of products for RT-PCR, qRT-PCR, microarray analysis, and reverse transcription is available separately from Invitrogen. For more information, visit our web site at www.invitrogen.com or contact Technical Support.

Troubleshooting

Observation	Cause	Solution
Low RNA yield	Incomplete lysis and homogenization	Use the appropriate method for lysate preparation based on your starting material as described on page 2. Decrease the amount of starting material used. Cut tissue samples into smaller pieces and ensure the tissue is completely immersed in the Lysis Buffer to achieve optimal lysis.
	Poor quality of starting material	The yield and quality of RNA isolated depends on the type and age of the starting material. Use fresh sample and process immediately after collection or freeze the sample at -80°C or in liquid nitrogen immediately after harvesting.
	Clogged RNA Spin Cartridge	Clear the homogenate and remove particulate or viscous material by centrifugation. Use only the supernatant for subsequent loading onto the Spin Cartridge.
	Ethanol not added to Wash Buffer II	Add ethanol to Wash Buffer II before use (page 2).
	Incorrect elution conditions	Add RNase-Free Water (30–100 μL) and perform incubation for 1 minute before centrifugation. To recover more RNA, be sure to use up to 3 sequential elutions of 100 μL each (3 \times 100 μL) Elution Buffer (refer to protocol on page 3).
RNA degraded	RNA contaminated with RNase	Follow the guidelines on page 2 to prevent RNase contamination.
Inhibition of downstream enzymatic reactions	Presence of ethanol in purified RNA	Traces of ethanol from the Wash Buffer II can inhibit downstream enzymatic reactions. Discard Wash Buffer II flow through. Place the Spin Cartridge into the Recovery Tube and centrifuge at $12,000 \times g$ for 1–2 minutes to completely dry the cartridge.
	Presence of salt in purified RNA	Use the correct order of Wash Buffers for washing. Always wash with Wash Buffer I followed by washing with Wash Buffer II.
Low A_{260}/A_{280} ratio	Sample was diluted in water	Use 10 mM Tris-HCl (pH 7.5) to dilute sample for OD measurements.

References

Chirgwin, J. M., Przybyla, A. E., MacDonald, R. J., and Rutter, W. Z. (1979) Isolation of Biologically Active Ribonucleic Acid from Sources Enriched in Ribonucleases. *Biochem. J.* 18, 5294-5299

Vogelstein, B., and Gillespie, D. (1979) Preparative and analytical purification of DNA from agarose. *Proc. Natl. Acad. Sci. USA* 76, 615-619

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008 or outlicensing@lifetech.com.

©2010 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.